

TattooStar R

Pictured on left

LASER TYPE

Q-Switched Ruby laser, class 4



The high power, optimized by the homogeneous spot, makes it possible to remove benign pigmented lesions quickly and effectively, without damaging the surrounding tissue. The pulses are absorbed selectively by the melanin pigments in the skin. These are pulverized by the laser energy and eliminated naturally by the lymphatic system or via the epidermis.



Pigmented lesions | before



Pigmented lesions | after

By courtesy of Dr. Said Hilton & Heike Heise, Dusseldorf, Germany

TattooStar removes a large variety of pigmented lesions:

